

Evaluation of the TGS TA system for the detection of anti-cytomegalovirus antibodies

Annalisa Cianflone, Maria Teresa Manco, Olivia Arpino, Alessia Paganini, Massimo De Paschale, Carlo Agrappi, Paola Mirri, Pierangelo Clerici

Microbiology Unit, Legnano Hospital, ASST Ovest Milanese, Legnano, Italy

Summary

Background and aims: The aim of the present study was to evaluate the new Technogenetics TGS TA system for detecting anti-Cytomegalovirus IgG and IgM antibodies and IgG avidity. The TGS TA system was compared with our routinely used system, LIAISON XL, for the detection of IgG and IgM antibodies. Only in positive IgM samples, TGS TA system was compared to an enzyme linked fluorescent assay (ELFA) test (VIDAS, BioMérieux, Marcy-l'Étoile, France) and with LIAISON XL system for the IgG avidity (if possible).

Materials and methods: Three hundred sera samples from pregnant women were examined with the TGS TA system and divided in 3 groups according to IgG and IgM screening LIAISON XL tests: 102 were non-immune women (Group 1), 98 were pregnant with past infection (Group 2) and 100 were pregnant with positive or equivocal IgM (95 with positive IgG and 5 with negative IgG) (Group 3).

Results: The overall concordance of the IgG results between LIAISON XL and TGS TA was 98.3%: 97.1% in Group 1, 100% in Group 2 and 98.0% in Group 3.

The overall concordance of the IgM results between LIAISON

XL and TGS TA was 92.1%: 100% in Group 1, 99.0% in Group 2 and 70.1% in Group 3. In Group 3, the concordance between the results of the IgG avidity with the LIAISON XL and TGS TA tests was 87.4%. Comparing the clinical diagnosis obtained with our protocol and that of the TGS TA system, the overall concordance was 94.3%: 97.1% in Group 1, 99.0% in Group 2 and 87.0% in Group 3.

Conclusions: In conclusion, the overall clinical concordance between the LIAISON XL/VIDAS protocol and the TGS TA system is excellent. TGA TA system shows to be a valuable tool able to clearly identify non-specific subjects, those with a non-recent infection and classify as either recent or past infection half of the subjects with undetermined infection with our protocol.

Introduction

The cytomegalovirus (CMV) is one of the leading agents responsible for congenital infections with clinical manifestations ranging from asymptomatic forms (majority of the cases) to serious damage to the foetus to, in rare cases, death by miscarriage (5,6,12,14,20). Screening campaigns have been suggested in pregnant, pre-pregnant women or newborn patients for preventive purposes in addition to health education campaigns but, in general, the scientific community is divided on the opportunity of carrying out serological screening campaigns (1,4,7, 8,13,15-18,21,22).

However, these tests are requested in any case by many general practitioners and gynaecologists and various analytical systems are currently available for this diagnosis (2,9-11,19) such as enzyme-linked immunosorbent assay (ELISA), chemiluminescence immunoassay (CLIA) and enzyme linked fluorescent assay (ELFA) which are able, depending on the instrument, to process a great number of samples in a short time. The commercial tests available generally have excellent sensitivity and specificity levels, with different formulations and the possible use of distinct antigens (natural, recombinant or synthetic one) which may be responsible for slight differences in the results obtained when comparing the available tests. The extent of these differences must be known in advance.

In this study we compared the CMV serology results obtained with our standard internal laboratory protocol, to the CMV results from the TGS TA system. In accordance with our laboratory protocol, on all samples with a positive or equivocal CMV IgM result obtained with LIAISON XL (DiaSorin, Saluggia, Italy), CMV IgM is also determined in a second round with VIDAS (BioMérieux, Marcy-l'Étoile, France) and a CMV IgG avidity analysis is performed on LIAISON XL.

Correspondence: Annalisa Cianflone Microbiology Unit, Legnano Hospital, ASST Ovest Milanese, via Papa Giovanni Paolo II, Legnano (MI), Italy.

Tel.: +39.0331.449319

E-mail: acianflone@libero.it

Key words: Cytomegalovirus; Pregnancy; IgM; IgG; Avidity.

Contributions: the authors contributed equally.

Conflict of interest: the authors declare no potential conflict of interest.

Received for publication: 18 January 2017.

Revision received: 15 March 2017.

Accepted for publication: 15 March 2017.

©Copyright A. Cianflone et al., 2017

Licensee PAGEPress, Italy

Microbiologia Medica 2017; 32:6581

doi:10.4081/mm.2017.6581

This article is distributed under the terms of the Creative Commons Attribution Noncommercial License (by-nc 4.0) which permits any non-commercial use, distribution, and reproduction in any medium, provided the original author(s) and source are credited.

Materials and Methods

Three hundred sera samples from women who came to our centre for the normal follow-up in pregnancy were examined and divided into the groups below:

Group 1: 102 samples from non-immune pregnant women

Group 2: 98 samples from pregnant women with past infection

Group 3: 100 samples from pregnant women with positive or equivocal IgM (95 with positive IgG and 5 with negative IgG).

Routinely, samples are tested using chemiluminescence LIAISON XL system for IgG and IgM (LIAISON CMV II IgG and CMV II IgM; DiaSorin, Italy). In positive subjects for IgM, samples are tested using an ELFA test (VIDAS Toxo IgM; BioMérieux, France). The avidity of the IgG (if possible) is determined with LIAISON XL (LIAISON CMV Avidity II; DiaSorin, Italy)

All samples are tested with TGS TA system for detecting IgG and IgM and, for Group 3, IgG avidity (TGS TA CMV IgG, TGS TA CMV IgM, TGS TA CMV IgG Avidity; Technogenetics, Milan, Italy).

For Groups 1 and 2, the TGS TA system was used on fresh samples immediately after the results with LIAISON XL test while, for Group 3, TGS TA was used on frozen samples stored on the basis of the results with our protocol.

The TGS TA CMV IgG and TGS TA CMV IgM kit employs an indirect two-step immunological method based on the principle of chemiluminescence. The solid phase is coated with a specific antigen and the conjugate is or an anti-human IgG antibody or anti-human IgG antibody labelled with an acridinium ester derivative. The TGS TA CMV IgG avidity uses a borate buffer able to prevalently break the binding Antigen (Ag)- Antibody (Ab) if the antibody is at low avidity.

The Interpretation scheme for the different tests are shown in Tables 1 and 2.

With reference to the IgG avidity, low avidity is strongly indicative of an infection occurring in the past 3 months for both systems; although high avidity does not exclude the possibility of a recent infection, it is strongly indicative of an infection at least 3 months before for both methods. The interpretation for moderate avidity is doubtful and the infection cannot be dated.

Results

The overall concordance and concordance in the three groups resulting from a comparison of the results obtained with the TGS TA and LIAISON XL systems is shown in Table 3.

Fifty-two samples were negative to the TGS IgM test in Group 3 (samples with positive or equivocal IgM for the LIAISON XL system). These samples were examined with the VIDAS IgM test and 23 were negative, 15 positive and 14 equivocal (Table 4).

In addition, the avidity was determined with both LIAISON XL and TGS TA in the 95 samples of Group 3 with IgG and the concordance was 87.4% (Table 5).

As far as the clinical assessment is concerned, the TGS TA system found 3 samples (2.9%) with low IgG titre (values 11.4-18.8 UA/mL) in Group 1 (non-immune women). All these subjects had values between 5.5 and 11.8 AU/mL for the LIAISON XL system. The clinical concordance is 97.1%.

The TGS TA test detected one case with positive IgM (not confirmed with VIDAS IgM) and low avidity with TGS TA and moderate with LIAISON XL in Group 2 (immune women). The clinical concordance is 99.0%.

Based on our protocol, patients in Group 3 (IgM positive or equivocal samples for LIAISON XL IgM) were classified as: 16 patients with a recent infection (samples confirmed with the VIDAS IgM and low avidity with LIAISON XL), of whom 3 without IgG; 14 with a non-determinable infection as they had moderate avidity with LIAISON XL (samples confirmed with the VIDAS IgM); 43 patients with past infection (samples confirmed with the VIDAS IgM and high avidity with LIAISON XL); 25 patients with IgM supposed non-specific not confirmed with the VIDAS IgM and high avidity with LIAISON XL; 2 non-immune patients, without IgG and IgM not confirmed with VIDAS IgM.

The interpretations obtained with the TGS TA system (IgG, IgM and IgG avidity) in these sub-groups are shown in Table 6. The concordance is 87.0%.

Considering all three groups, the clinical concordance is shown in Table 7.

Table 1. Interpretation scheme for IgG and IgM anti-cytomegalovirus with LIAISON XL systems, VIDAS and TGS TA.

Interpretation scheme	IgG		IgM		
	LIAISON XL, IU/mL	TGS TA, AU/mL	LIAISON XL, AU/mL	VIDAS (index)	TGS TA, AU/mL
Negative	<12	<10	<18	<0.70	<6
Equivocal	12-14	-	18-22	0.70-0.90	-
Positive	>14	≥10	>22	>0.90	≥6

Table 2. Interpretation scheme for the anti-cytomegalovirus IgG avidity tests with the LIAISON XL and TGS TA.

Interpretation scheme	IgG avidity	
	LIAISON XL (index)	TGS TA (index)
Low	<0.150	<0.150
Moderate	0.150-0.250	0.150-0.200
High	>0.250	≥0.200

Discussion and Conclusions

The concordance between LIAISON XL and TGS TA in the IgG tests is excellent in the three groups. The few discrepancies are due to the IgG values close the cut-off. The concordance for the IgM tests is excellent in the first two groups; in the third group, the discordance is greater, also because there were subjects with non-specific IgM in this group. Subjects with IgM positive in LIAISON XL (also with VIDAS) and negative in the TGS TA test, the avidity, in all cases except one, was high with both LIAISON XL and TGS TA. It should be recalled that the TGS TA system for IgM does not have a grey doubtful area around the cut-off, unlike the

LIAISON XL and VIDAS, therefore the calculation of the concordance suffers from this limitation.

The concordance between LIAISON XL and TGS TA in the IgG avidity tests is good. Most of samples with moderate avidity with LIAISON XL test, shown a high or low avidity values with TGS TA system. Therefore, the indeterminate infections found with our protocol could be defined as recent or past with TGS TA method.

A limit of the study design is the lack of cases of recent CMV infection due to low incidence of acute infection in pregnant women in our area (3) that makes quite rare the observation of recent infection rare. Another limitation of this study is the selection of samples in Group 3 using our protocol which provides screening for IgG and IgM with LIAISON XL. On the other hand,

Table 3. Overall concordance and concordance in the three groups between LIAISON XL and TGS TA for the detection of anti-cytomegalovirus IgG and IgM in samples from pregnant women.

	IgG		IgM	
	TGS TA neg	TGS TA pos	TGS TA neg	TGS TA pos
Overall concordance in 300 samples from pregnant women*				
LIAISON XL neg	102	5	199	1
LIAISON XL pos	0	193	20	47
LIAISON XL eq	-	-	32	1
Group 1: concordance in 102 samples from non-immune pregnant women**				
LIAISON XL neg	99	3	102	0
LIAISON XL pos	0	0	0	0
Group 2: Concordance in 98 samples from pregnant women with past infection***				
LIAISON XL neg	0	0	97	1
LIAISON XL pos	0	98	0	0
Group 3: concordance in 100 samples from pregnant women with positive or equivocal IgM with LIAISON XL****				
LIAISON XL neg	3	2	0	0
LIAISON XL pos	0	95	20	47
LIAISON XL eq	-	-	32	1

*Overall concordance (excluding equivocal cases): IgG 98.3%; IgM 92.1%. **Overall concordance: IgG 97.1%; IgM 100%. ***Overall concordance: IgG 100%; IgM 99%. ****Overall concordance (excluding equivocal cases): IgG 98%; IgM 70.1%.

Table 4. Results of the VIDAS on 52 samples equivocal or positive to the LIAISON XL anti-cytomegalovirus IgM and negative with the TGS TA IgM.

IgM positive or equivocal for the LIAISON XL system	LIAISON XL	VIDAS	TGS TA
Negative	0	23	52
Positive	20	15	0
Equivocal	32	14	0

Table 5. Concordance between the LIAISON XL and TGS TA anti-cytomegalovirus IgG avidity tests in 95 individuals positive or equivocal (with IgG positive) for LIAISON XL IgM.

Avidity	TGSTGS TA avidity			Total
	Low	Moderate	High	
LIAISON XL				
Low	11	1	1	13
Moderate	5	5	4	14
High	0	1	67	68
Total	16	7	72	95

Concordance: 87.4%.

Table 6. Comparison between the clinical interpretations given with cytomegalovirus LIAISON XL/VIDAS and TGS TA systems.

LIAISON XL/VIDAS system	Probable recent infection	Undeterminate infection	TGS TA system Probable non-recent infection	Non immune	Total
Probable recent infection	13	1	1	1	16
Undeterminate infection	5	5	4	0	14
Probable non-recent infection	0	1	67	0	68
Non immune	0	0	0	2	2
Total	18	7	72	3	100

Concordance 87.0%.

Table 7. Overall clinical agreement between the results obtained with the cytomegalovirus LIAISON XL/VIDAS protocol and those of the TGS TA system.

LIAISON XL/VIDAS protocol	Probable recent infection	Undeterminate infection	TGS TA system Probable non-recent infection	Non immune	Total
Probable recent infection	13	1	1	1	16
Undeterminate infection	6	5	4	0	15
Probable non-recent infection	0	1	164	0	165
Non immune	0	0	3	101	104
TOTAL	19	7	172	102	300

Total concordance 94.3%.

our protocol leads to the definition of the type of infection that was compared with the interpretation given by the TGS TA system. The clinical interpretation correlates well between our protocol and TGS TA system. Taking into account all three groups, the TGS TA system clearly identifies many subjects non-specific and those with non-recent infection. Further, although it missed two cases of recent infection in Group 3, it found one more in Group 2 compared to our protocol. In addition, more than half of samples with undetermined infection with our protocol are clearly defined by TGS TA as either recent or past infection.

In conclusion the overall clinical concordance between the two systems is good and the TGS TA is a useful system to evaluate the immune status related to CMV infection.

References

- Ahlfors K, Ivarsson SA, Harris S. Report on a long-term study of maternal and congenital cytomegalovirus infection in Sweden. Review of prospective studies available in the literature. *Scand J Infect Dis* 1999;31:443-57.
- De Paschale M, Agrappi C, Manco MT, Clerici P. Positive predictive value of anti-HCMV IgM as an index of primary infection. *J Virol Methods* 2010;168:121-5.
- De Paschale M, Agrappi C, Manco MT, et al. Incidence and risk of cytomegalovirus infection during pregnancy in an urban area of Northern Italy. *Infect Dis Obstet Gynecol* 2009;206505.
- Enders G, Bäder U, Lindemann L, et al. Prenatal diagnosis of congenital cytomegalovirus infection in 189 pregnancies with known outcome. *Prenat Diagn* 2001;21:362-77.
- Fowler KB, Stagno S, Pass RF, et al. The outcome of congenital cytomegalovirus infection in relation to maternal antibody status. *N Engl J Med* 1992;326:663-7.
- Fowler KB, Stagno S, Pass RF. Maternal immunity and prevention of congenital cytomegalovirus infection. *JAMA* 2003;289:1008-11.
- Gratacap-Cavallier B, Morand P, Dutertre N, et al. Cytomegalovirus infection in pregnant women. Seroepidemiological prospective study in 1,018 women in Isere. *J Gynecol Obstet Biol Reprod (Paris)* 1998;27:161-6.
- Hagay ZJ, Biran G, Ornoy A, Reece EA. Congenital cytomegalovirus infection: a long-standing problem still seeking a solution. *Am J Obstet Gynecol* 1996;174:241-5.
- Kaneko M, Ohhashi M, Minematsu T, et al. Maternal immunoglobulin G avidity as a diagnostic tool to identify pregnant women at risk of congenital cytomegalovirus infection. *Infect Chemother* 2017;23:173-6.
- Lazzarotto T, Gabrielli L, Lanari M, et al. Congenital cytomegalovirus infection: recent advances in the diagnosis of maternal infection. *Hum Immunol* 2004;65:410-5.
- Lazzarotto T, Guerra B, Gabrielli L, et al. Update on the prevention, diagnosis and management of cytomegalovirus infection during pregnancy. *Clin Microbiol Infect* 2011;17:1285-93.
- Leruez-Ville M, Ville Y. Fetal cytomegalovirus infection. *Best Pract Res Clin Obstet Gynaecol* 2017;38:97-107.
- Mazeron MC. Conséquence néonatales des maladies sexuellement transmises. Conduite à tenir devant une infection à cytomegalovirus. *Med Mal Infect* 1994;24:485-95.
- Mestas E. Congenital cytomegalovirus. *Adv Neonatal Care* 2016;16:60-5.
- Peckham CS, Chin KS, Coleman JC, et al. Cytomegalovirus infection in pregnancy: preliminary findings from a prospective study. *Lancet* 1983;i:1352-3.
- Revello MG, Gerna G. Diagnosis and management of human cytomegalovirus infection in the mother, fetus, and newborn infant. *Clin Microbiol Rev* 2002;15:680-715.
- Revello MG, Tibaldi C, Masuelli G, et al. Prevention of pri-

- mary cytomegalovirus infection in pregnancy. *EBioMedicine* 2015;2:1205-10.
18. Schoub BD, Johnson S, McAnerney JM, et al. Is antenatal screening for rubella and cytomegalovirus justified? *S Afr Med J* 1993;83:108-10.
 19. Sellier Y, Guilleminot T, Ville Y, et al. Comparison of the LIAISON(®) CMV IgG Avidity II and the VIDAS(®) CMVIgG Avidity II assays for the diagnosis of primary infection in pregnant women. *J Clin Virol* 2015;72:46-8.
 20. Stagno S, Pass RF, Cloud G, et al. Primary cytomegalovirus infection in pregnancy : Incidence, transmission to fetus and clinical outcome. *JAMA* 1986;256:1904-8.
 21. Thackeray R, Magnusson BM. Women's attitudes toward practicing cytomegalovirus prevention behaviors. *Prev Med Rep* 2016;4:517-24.
 22. Yow MD, Demmler GJ. Congenital cytomegalovirus disease – 20 years is long enough. *N Engl J Med* 1992;326:702-3.

Non-commercial use only